

IDENTIFY  *chrome*  
Evolution of Light



## The Chrome Advantage

The Brightest, Lightest Headlight for the Medical Profession

Brought to you by *micronix*

(877) 379-7340 • [www.MicronixSystems.com](http://www.MicronixSystems.com)

## Intelligent Design

Identify headlights bring the latest LED illumination technology to your medical procedures. Emitting a cool, yet powerful white light, they more than double the brightness of standard LEDs or conventional headlights. At **250,000 LUX**, Identify Chrome provides maximum, low-temperature operating brightness without the undesirable heat generated by traditional halogen headlights. The Identify Chrome shines a beam from 9cm – 20cm in size, with the broadest range of superior illumination.

Identify Chrome features maximum adjustability in the controls, including digital brightness, dual axle beam angle and spot diameter for ultimate performance.

## Cordless Advantage

Powerful streamlined battery technologies integrated into the headstrap enable total freedom of movement. This eliminates the need for bulky belt packs or restrictive tethering cables. The six lithium polymer batteries provide 12 hours of wireless power at full brightness and can be hot-swapped for infinite periods of use.

## Maximum Comfort

Weighing just 11 ounces, Identify Chrome is one of the lightest headlights available for the medical profession. The Identify Chrome Slimstrap™ minimizes fatigue and delivers all day, ultra-comfortable wear.

Identify Chrome interchangeable bezels for specific applications and a detachable camera accessory for recording procedures are also available.



**Each Hospital Set Includes a Headlight, 6 Lithium Polymer Battery Packs, 4-Bay Desktop Charger, a USB Cord, Power Outlet Adapter, and a Secure Carrying Case.**



Manufactured by RIVERPOINT MEDICAL  
 ISO 13485:2003;  
 ISO 9001:2008;  
 FDA Registered



## SPECIFICATIONS

### Headlight Specification



<b>BRIGHTNESS</b>	Up to 250,000 LUX with MC7 Bezel	
<b>SPOT SIZE</b>	9 cm – 20 cm	●●
<b>DIGITAL DIMMING CONTROL</b>	Yes	✓
<b>BATTERIES</b>	6 Lithium Polymer	000000
<b>RUN TIME</b>	12 Hours with 6 Batteries	○
<b>HOT SWAPPING POSSIBLE</b>	Yes	✓
<b>CHARGE</b>	Micro USB Port – Charge from Computer or Power Outlet	
<b>4-BAY DESKTOP CHARGER</b>	Included	✓
<b>BATTERY GAUGE</b>	Yes	✓
<b>ERGONOMICS</b>	Ultra-Comfort Chrome Slimstrap™ -Integrated Batteries, Digital Adjustment	
<b>WEIGHT</b>	11 ounces	
<b>CAMERA OPTION</b>	Yes	✓
<b>INTERCHANGEABLE BEZELS</b>	Available	✓

### Rechargeable Batteries (6 included in set)

<b>TYPE</b>	Lithium Polymer
<b>VOLTAGE</b>	3.7V
<b>RUN TIME</b>	12 Hours with 6 Batteries (at Full Brightness) 5 Indicator LEDs Display Battery Levels Infinite Run Time with Hot swap

**Ask for your free trial of this sleek wearable lighting solution for the operating room.**

**(877) 379-7340**

Call for more Information or visit our website at:

[www.MicronixSystems.com](http://www.MicronixSystems.com)

44 Commerce Street  
 Springfield, New Jersey 07081  
[Sales@MicronixSystems.com](mailto:Sales@MicronixSystems.com)