

IDENTIFY  *surgical*
Evolution of Light



The Surgical Advantage

The Brightest, Lightest Headlights for the Medical Profession

Brought to you by *micronix*

Intelligent Design

Identify headlights bring the latest LED illumination technology to your medical procedures. Emitting a cool, yet powerful white light, they more than double the brightness of standard LEDs or conventional headlights. At **200,000 LUX**, Identify Surgical provides superior, low-temperature operating brightness without the undesirable heat generated by traditional halogen headlights. The Identify Surgical shines a beam from 3.5 cm – 8.5 cm, an ideal size for strong light with optimum depth illumination.

Identify Surgical features maximum adjustability in the controls, including digital brightness, dual axle beam angle and diameter for ultimate performance.

Cordless Advantage

Powerful streamlined battery technologies integrated into the headstrap enable total freedom of movement. This eliminates the need for bulky belt packs or restrictive tethering cables. The two lithium polymer batteries provide six to eight hours of wireless power at full brightness and with additional batteries and charger, can be hot-swapped for infinite periods of use.

Maximum Comfort

Weighing just 11 ounces, Identify Surgical is one of the lightest headlights available for the medical profession. Ergonomically balanced, Identify Surgical has a cushioned and easily adjustable headband that maximizes wearing comfort and minimizes fatigue.



Each Standard Set Includes a Headlight, 2 Lithium Polymer Battery Packs, a USB Cord, Power Outlet Adapter, and a Secure Carrying Case.



Manufactured by RIVERPOINT MEDICAL
 ISO 13485:2003;
 ISO 9001:2008;
 FDA Registered



SPECIFICATIONS

Headlight Specification



BRIGHTNESS	Up to 200,000 LUX	
SPOT SIZE	3.5 cm – 8.5 cm	• ●
DIGITAL DIMMING CONTROL	Yes	✓
BATTERIES	2 Lithium Polymer	🔋
RUN TIME	6-8 Hours with 2 Batteries	🕒
HOT SWAPPING POSSIBLE	Yes	✓
CHARGE	Micro USB Port – Charge from Computer or Power Outlet	
4-BAY DESKTOP CHARGER	Available	
BATTERY GAUGE	Yes	✓
ERGONOMICS	Comfortable, Cushioned Headband-Integrated Digital Adjustment	
WEIGHT	11 ounces	
CAMERA OPTION	No	
INTERCHANGEABLE BEZELS	No	

Rechargeable Batteries (2 included in set)

TYPE	Lithium Polymer
VOLTAGE	3.7V
RUN TIME	6+ Hours with 2 Batteries (at Full Brightness) 5 Indicator LEDs Display Battery Levels Infinite Run Time with Hot Swap

Ask for your free trial of this sleek wearable lighting solution for the operating room.

(877) 379-7340

Call for more Information
 or visit our website at:

www.MicronixSystems.com

44 Commerce Street
 Springfield, New Jersey 07081
Sales@MicronixSystems.com